# 4655\_CEPILLO DIENTES BESOL

## FRONT PART

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## BACK PART

### 4655\_Toothbrush Besol User manual (Eng)

It is recommended to read this user manual carefully and keep it for future reference.

#### Operation:

- This electric toothbrush has a head specially designed to penetrate into the inter-dental spaces while respecting gums. The flexibility of its bristle allows an effective teeth brushing.
- To avoid splatters, place the brush on your teeth before starting up the toothbrush. Then guide the brush slowly from a tooth to another one, insting a few seconds on the surface of every tooth. Brush the outside, then the inside of teeth and finally, the chewing surfaces. Do not press too hardly, do not rub, let simply the toothbrush does the work by itself.
- The brushing time recommended by the dentist is about 3 minutes.

#### Recommendation:

- The cleaning and maintenance must not be carried out by children without adult supervision.
- Children must not play with this device.
- Apart from the batteries, do not replace any element of this item; In the event of any piece being damaged (such us the casing or the head), it is advisable to dispose of the item.

- Do not immerse the body of the toothbrush in water or any other liquid.
- This product operates with 2 AA/LR6 batteries.
- Insert the batteries according to the polarity signs indicated.
- The batteries not rechargeable must not be recharged.
- Do not cut or damage the battery ends.
- Do not mix new batteries with use dones or different types of batteries and different brands.
- Remove exhausted batteries and dispose of them appropriately.
- Remove batteries before disposing the device.
- Remove batteries when item is not going to be used for a long period of time.

This symbol both in the product and its user manual indicates that at the end of the devices electrical life, it should be recycled separately from your household waste. There are adequate means to dispose of this material for proper recycling. For more information, contact your local authorities.

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